

G6: Focused on the patient experience



- Size Reduction
- 1x per Day Calibration
- 10 Day Use Life
- Interference Blocking
- Performance Improvement

1) Represents calibrations required

2) Mean Absolute Relative Difference (MARD) between CGM readings and blood glucose readings

3) Proportion of sensor values within ± 20 mg/dL of the YSI reference value for glucose levels ≤ 80 mg/dL and $\pm 20\%$ at YSI glucose levels > 80 mg/dL

* Note: G6 data is pre-pivotal

		G5 2/day¹	G6 1/day¹ ALL DATA	G6 1/day¹ ADULTS	G6 No Cal¹ ALL DATA
Overall	N	2263	993	646	1009
	MARD ²	9%	8.1%	7.2%	8.8%
	%20/20 ³	93%	96.1%	99.4%	95.5%
Day 1	N	680	161	95	165
	MARD ²	10.7%	10.4%	7.3%	9.8%
	%20/20 ³	84%	90.1%	100%	94.5%